3M

Hygiene mask

Maschera igienica

Mascarilla de higiene

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Mundschutz

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2

3

Fitting Instructions: 1. Cup Mask under chin with nose clip up. 2. Pull headband up to crown of head. 3. Press the nose clip around the nose.

Warnings: I) This mask is **NOT A RESPIRATORY PROTECTIVE DEVICE** and should not be used to protect your health. II) It should not be used in **ANY** potentially hazardous atmospheres. For example atmospheres containing: dusts, mists, fibres, gases or vapours, solvents or paintspray may be hazardous to health. III) Leave the area if breathing becomes difficult, dizziness or other distress occurs. IV) Never alter or modify this device. V) Use only in well ventilated atmospheres containing sufficient oxygen to support life. VI) Misuse may result in sickness or permanent disability.

Information & Instructions: This mask **SHOULD ONLY** be used to assist in preventing pollution of the environment or the process by the exhaled breath of the wearer. **Hygiene masks are not defined as Personal Protective Equipment** as described in European Directive 89/686 and **SHOULD NOT** be used where a device complying with European Directive 89/656 is required. **Hygiene masks** (also called Comfort Masks) **are defined as any mask** (or gauze pad) from any supplier which **is not designed to meet the requirements of the relevant European Standard** (e.g. EN149). 3M has not recommended the use of Hygiene masks as Personal Protective Equipment. The above statement does not apply to certified disposable respirators that do meet relevant standards. Phone the 3M Health & Safety Helpline on 0870 6080060 (UK) or 1 800 320500 (Ireland) for advice on selection of 3M's range of certified respirators.

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